



DISTRIBUTED BY



1-800-361-5361

Alfa Laval LKH Prime UltraPure

Centrifugal pumps

Introduction

The Alfa Laval LKH Prime UltraPure Centrifugal Pump is designed for use in high-purity applications where high efficiency, exceptional cleanability, contamination safety, robust design and low maintenance are of paramount importance.

Precision-engineered, the LKH Prime UltraPure delivers greater energy efficiency than similar pumps. Its optimized design, premium motor, tight tolerances and advanced impeller and airscrew design minimize recirculation and reduce energy consumption.

Applications

The Alfa Laval LKH Prime UltraPure is designed to meet the stringent demands and regulations of high-purity applications across the biotechnology and pharmaceutical industries that require equipment with the highest material integrity. It is ideal for tank emptying and CIP return applications; it has verified and effective CIP cleanability. The LKH Prime UltraPure can also be used as a product pump.

All pumps are delivered with a complete Alfa Laval Q-doc package. Q-doc provides easier validation, proof of origin and compliance for inspection purposes according to Good Manufacturing Practice (GMP) and ASME BPE requirements.

The LKH Prime UltraPure pump is available in two sizes to handle capacities up to 70 m³/h and differential pressures up to 4 bar at 50 Hz.



Benefits

- Energy efficient: superior efficiency resulting in reduced energy consumption and CO₂ footprint.
- Quiet: operates very quietly compared to other self-priming pumps, thereby improving the working environment.
- Low contamination risk: comes with full material traceability and USP Class VI elastomers to reduce risk of process contamination from extractables.
- Smooth qualification, validation and process control: material traceability, and pump supplied with the Alfa Laval Q-doc package in line with Good Documentation Practices (GDP).

Standard design

All media contacting steel components like pump casing, impeller, airscrew, front cover, recirculation pipe and backplate are in W. 1.4404 (AISI 316L) with material traceability 3.1 according to EN 10204. Product wetted elastomers are specified to USP Class VI, 249.8°F, Chapter 88 and Chapter

87. A stainless steel shroud protects the motor and four adjustable stainless steel legs support the complete unit.

A compression coupling securely attaches the stub shaft to the motor shaft with precision alignment, and the semi-open impeller with a special vane design ensures efficient handling of the product as it moves through the pump.

As standard, the LKH prime pump is equipped with a single mechanical shaft seal but is also available with a double mechanical shaft seal. The front-loading shaft seal, with the spring and washers mounted on the atmospheric side, makes maintenance fast, easy and inexpensive. It takes just a few minutes to replace the shaft seal. In addition, the balanced design minimizes the risk of seal opening during unforeseen pressure shock.

Working principle

On applications where the pumped media contains a mixture of air and liquid in the suction line, airscrew rotation causes the formation of a continuous liquid ring within the canister.

Due to the eccentric position of the canister relative to the airscrew, an air chamber forms between the liquid ring and the airscrew, which separates into air pockets between the air-screw vanes.

The continuous rotation of the airscrew forces air pockets through the canister into the suction stage of the impeller which are then pumped out via the discharge.

Liquid is returned from the discharge via the recirculation pipe into the canister to ensure the liquid ring is maintained at all times. When there is no air present, the canister and recirculation loop have no function and are fully filled with liquid. The liquid passes through the canister into the suction stage of the impeller, allowing the pump to act as a traditional centrifugal pump.

Certificates



Authorized to carry the 3A symbol

TECHNICAL DATA

Standard materials	
Product wetted steel parts:	AISI 316L and 329L with material traceability 3.1 acc. to EN 10204 (Mill test reports).
Other steel parts:	Stainless steel.
Inside surface finish:	EP Ra ≤ 15 µin.
Product wetted elastomers:	EPDM- USP Class VI, 249.8°F Chapter 88, and Chapter 87.

Motor

Standard C-faced, foot mounted motor according to NEMA standard. 3500 RPM. Premium efficiency, Class F. Note different frame sizes.

Connections

Connections for double mechanical shaft seal:	1/8" NPT.
---	-----------

Min/max motor speed

Air evacuation:	2800 - 3600 rpm.
Pumping product (no air):	900 - 3600 rpm.


OPERATING DATA

Pressure	
Max inlet pressure:	72.5 PSI (5 bar).

Temperature	
Temperature range:	14°F to 284°F (EPDM).

Double mechanical shaft seal	
Water pressure inlet:	Max. 72.5 PSI (5 bar).
Water consumption:	4-8 US gph.

Dimensions



Pump specific measures

Pump Model	LKH Prime UltraPure 10	LKH Prime UltraPure 20
A	6.85	7.36
B	3.35	3.62
C	8.74	9.76
D	9.72	9.96
E	9.64	11.02
P1	4.57	4.84
P2	3.24	3.27
T1	8.63	9.41
T2	6.35	6.69

Motor specific measures

Motor TC/TSC	182TC	184TC	213TC	215TC	254TC	256TC
Motor HP	3.0	5.0	7.5	10.0	15.0	20.0
F(max) ¹	9.09	9.09	9.84	9.84	10.87	10.87
G	4.49	4.49	5.39	5.39	6.77	6.77
H	11.10	11.10	13.07	13.07	17.80	17.80
I	17.99	18.15	21.02	21.02	26.06	27.83

¹ Possible to reduce dimension F by min. 2.32 inch for all pump models.

Motor overview

Pump Model	LKH Prime UltraPure 10	LKH Prime UltraPure 20
Motor range (TC/TSC)	182TC-215TC	182TC-256TC



Note! Dimensional data are based on 2 pole, Sterling motors.

Connections

Pump Model		LKH Prime UltraPure 10	LKH Prime UltraPure 20
TRI-Clamp	M1	1.13	1.13
	M2	1.13	1.13
J1 ¹		2.00"	2.50"
J2 ¹		2.00"	2.00"

¹ Other dimensions available on request.

	TC Clamp
1/2"	0.5
3/4"	0.75

Flow chart




Figure 1. Frequency: 60Hz - Speed (synchr): 3600 rpm

Options

- Impeller with reduced diameter.
- Motor with increased safety/flame proof motor.
- Double mechanical shaft seal.
- Product wetted surface finish mechanically polished to Ra 20 μ m.
- Passivated surface.
- Product wetted elastomers FPM or FEP to USP Class VI, 249.8°F Chapter 88, and Chapter 87.
- 3/4" drain connection.
- No drain.
- 0° outlet position.
- Hydrostatic testing with certificate.
- Surface finish measurement with certificate.

Q-doc

Standard documentation package:

- Declaration of compliance to EN 10204 type 3.1 (MTR).
- Declaration of compliance to the U.S. Food & Drug Administration CFR 21 (non-metallic parts).
- Declaration of compliance to the U.S. Pharmacopeia (Elastomers and polymers).
- TSE (Transmissible Spongiform Encephalopathy) / ADI (Animal Derivative Ingredient) declaration.
- Declaration of surface finish compliance.
- Declaration of passivation and electro polishing (if specified).
- 3.1 certification in accordance to EN10204.
- Pump performance test certificate.

Optional documentation:

- Hydrostatic test certificate.
- Surface measurement report.

Ordering

Please state the following when ordering:

- Pump size.
- Connections.
- Impeller diameter.
- Motor size.
- Voltage and frequency.
- Flow, pressure and temperature.
- Density and viscosity of the product.
- Options.

This document and its contents are subject to copyrights and other intellectual property rights owned by Alfa Laval Corporate AB. No part of this document may be copied, re-produced or transmitted in any form or by any means, or for any purpose, without Alfa Laval Corporate AB's prior express written permission. Information and services provided in this document are made as a benefit and service to the user, and no representations or warranties are made about the accuracy or suitability of this information and these services for any purpose. All rights are reserved.

200006080-1-EN-US

© Alfa Laval Corporate AB

How to contact Alfa Laval

Up-to-date Alfa Laval contact details for all countries are always available on our website at www.alfalaval.com

DISTRIBUTED BY
HARCO
ENTERPRISES LTD.
1-800-361-5361

